

The seamless knitted
tubular bandage
for universal use:



Stülpa®

Stülpa is a seamless knitted tubular bandage with high two-way stretch. Not only does it provide very secure dressing retention, but the material is also particularly ideal for lining and covering zinc paste and plaster of Paris bandages, as well as for immobilisation or extension bandages.



Stülpa widens when stretched width-wise and narrows when stretched length-wise. This interaction is the basis of the benefits of the Stülpa bandaging technique. Stülpa is pulled over the part of the body to be bandaged and anchored appropriately according to the type of dressing. This provides a permanently non-slip finish.

Description

▶ Their high elasticity makes Stülpa bandages quick and easy to apply without additional fixation.

▶ They fit all body parts without creasing, slipping or causing a tourniquet effect. Joint movements are not impaired.

▶ Unlike other finished bandages, Stülpa bandages provide a smooth surface. This means that they do not loosen even under mechanical stress. Wounds or sensitive parts of the skin are reliably protected.

▶ Stülpa is made of 70% viscose and 30% cotton, making it pleasant to wear, absorbent, and air-permeable; it may also be sterilised.

▶ Stülpa can be cut anywhere without fraying or mesh run.

▶ Stülpa rolls are available in nine sizes with tube widths ranging from 1.5 to 24 cm. Ready-for-use Stülpa bandages are also available in the most commonly used sizes.

Uses

For all types of retention bandages, immobilisation bandages (Desault's bandages) and extension bandages, for lining and covering zinc paste or plaster of Paris bandages, cushioned dressings, for covering splints and under cast padding.

Presentations

		Unit	Code No.	Units per Case
15 m rolls	Size 0 R	1 roll	427 110	20
	Size 1 R	1 roll	427 111	24
	Size 2 R	1 roll	427 112	32
	Size 3 R	1 roll	427 113	28
	Size 4 R	1 roll	427 114	14
	Size 5 R	1 roll	427 115	12
6 m rolls	Size 6 R	1 roll	427 146	10
	Size 7 R	1 roll	427 147	16
	Size 8 R	1 roll	427 148	10

Ready-for-use bandages

individually wrapped	Size 1	1 piece	427 211	20/240
	Size 3	1 piece	427 213	10/240
	Size 4	1 piece	427 214	10/240
in cardboard boxes	Size 1	50 pieces	427 261	12
	Size 3	10 pieces	427 253	33
	Size 4	10 pieces	427 254	24

Size table

- Size 0 = Children's finger and toe bandages
- Size 1 = Finger bandages
- Size 2 = Arm and children's leg bandages
- Size 3 = Foot, leg, children's head and axilla bandages
- Size 4 = Head, leg and axilla bandages, face masks
- Size 5 = Children's trunk, head, thigh and axilla bandages, extension bandages
- Size 6 = Children's trunk and thigh bandages
- Size 7 = Desault's, trunk and extra-large thigh bandages
- Size 8 = Desault's and trunk bandages