



# The powder-free examination gloves at a glance.

	Peha-soft® pf	Peha-soft® satin pf	Peha-soft® syntex pf	Peha-soft® vinyl pf	Peha-fol® pf	Digitil® N pf
Latex	●	●				
Latex-free			●	●	●	●
Sterile	●					
Non-sterile	●	●	●	●	●	●
Textured fingertips	●					●
Smooth surface		●	●	●	●	



# Guidelines for glove use.

	Gloves made of natural latex		
	Peha-soft® pf	Peha-soft® pf sterile	Peha-soft® satin pf
<b>Applications in the operating theatre area:</b>			
Examination procedures	●		●
Sterilization			●
Anaesthesia	●		●
Circulating nursing activities	●		●
Skin disinfection		●	
<b>Applications in hospitals/physicians practices/nursing care units:</b>			
Handling patients	●	●	●
Handling contaminated material	●		●
Basic medical care			●

	Gloves for persons allergic to natural latex			
	Peha-soft® syntex pf	Peha-soft® vinyl pf	Peha-fol® pf	Digitil® N pf
	●			●
				●
	●	●		●
	●	●		●
	●			●
				●
		●	●	

If a glove suitable for multiple applications, the choice of examination glove should be made dependent on the actual situation.